

# Ready-to-use Human Arrays

*at your fingertips*

## UHN Human Microarrays



- ✓ Ready-to-use high-density spotted arrays
- ✓ “Face-to-face” hybridization of any 2 expression arrays
- ✓ Reliable and Affordable
- ✓ All arrays are available through our Expression Service



## Substantial Coverage of the Human Genome

- 19K - most popular set used in gene expression and profiling studies
- 8K - well-characterized clones with OMIM, RefSeq or SwissProt annotations
- 27K - combined set of 19K and 8K
- 10K - combined set of 8K and the retired 1.7K set (a subset of 19K)
- 12K - CpG-island clones for gene regulation studies (The Sanger Institute)
- 34.6K - 70-mer oligonucleotides based on Ensembl database

Human Clone Collection <sup>1</sup>					
Clone Sets	H12K CpG	H1.7K	H8K	H19K	oH34.6K
No. of Clones	12,192	1,718	8,196	19,008	34,580
Array Types	HCGI12K	H1.7K	H8K	H19K	oH34K
		H10K			



1: Annotations subject to change; data based on Unigene database (Oct 24/2003)

2: 509 clones contain accession annotations only

UHN Microarray Centre  
 University Health Network  
 200 Elizabeth Street  
 Toronto, Ontario M5G 2C4  
 Tel: 416-340-4259 Fax: 416-340-4004  
 Toll: 877-294-4410 (Canada & USA only)  
**Visit us at [www.microarrays.ca](http://www.microarrays.ca)**

